

CHIRURGISCHE FEDERÖHRNADELN

CHIRURGISCHE FEDERÖHRNADELN im praktischen Schiebeblistert:

Durch die Einzelanordnung ist keine gegenseitige Beschädigung der Nadelspitzen mehr möglich

- Schnelle und praktische Einzelentnahme aus dem unsterilen Blister
- Sichere und saubere Lagerung durch einfachen Wiederverschluß

Gehärteter Edelstahl	■ hohe Biegefestigkeit
	■ garantiert rostfrei
Optimiertes Federöhrdesign	■ sicherer Fadenhalt
	■ vermindert das Risiko einer Verletzung der Fadenoberfläche
	■ schnelles und leichtes Einfädeln durch optimale Elastizität des Federöhres
Präzisionsschliff	■ Reduzierung des Wundtraumas durch scharfen Kantenschliff
	■ hohe Penetrationsfähigkeit

SPRING-EYED NEEDLES

SPRING-EYED surgical NEEDLES in the practical sliding blister pack











Mutual damage to the needle points no longer possible due to the individual packing

- *Fast and practical removal of individual needles*
- *Safe and clean storage in simple reclosable packs*

<i>Tempered stainless steel</i>	■ <i>high bending resistance</i>
	■ <i>guaranteed rustproof</i>
<i>Optimised spring eye design</i>	■ <i>secure hold for the suture</i>
	■ <i>reduces the risk of damage to the suture's surface</i>
	■ <i>fast and easy suture insertion due to the optimal elasticity of the spring eye</i>
<i>Precision honing</i>	■ <i>sharply honed cutting edge reduces tissue trauma,</i>
	■ <i>excellent penetration ability</i>








CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL

SPRING-EYED NEEDLES – NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
5/8-Kreis Rundkörper, standard = H(318) <i>5/8 circle, round bodied, standard = H(318)</i>				
● FR 42 	H 318 / 2	FR 42 (130)	OE-FR 42	3
1/2-Kreis Rundkörper, standard = G(312) <i>1/2 circle, round bodied, standard = G(312)</i>				
● HR 20 	G 312 / 15	HR 20(70)	OE-HR 20	6
● HR 23 	G 312 / 14	HR 23(80)	OE-HR 23	6
● HR 26 	G 312 / 13	HR 26(80)	OE-HR 26	6
● HR 30 	G 312 / 12	HR 30(80)	OE-HR 30	6
● HR 34 	G 312 / 11	HR 34(90)	OE-HR 34	6
● HR 38 	G 312 / 10	HR 38(90)	OE-HR 38	6
● HR 42 	G 312 / 9	HR 42(100)	OE-HR 42	6
● HR 46 	G 312 / 8	HR 46(100)	OE-HR 46	6
● HR 50 	G 312 / 7	HR 50(110)	OE-HR 50	6











CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL

SPRING-EYED NEEDLES – NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
<p>● HR 54</p> 	G 312 /6	HR 54(110)	OE-HR 54	3
<p>● HR 58</p> 	G 312 /5	HR 58(120)	OE-HR 58	3
<p>● HR 63</p> 	G 312 /4	HR 63(130)	OE-HR 63	3
<p>● HR 68</p> 	G 312 /3	HR 68(130)	OE-HR 68	3
<p>● HR 73</p> 	G 312 /2	HR 73(140)	OE-HR 73	3
<p>● HR 80</p> 	G 312 /1	HR 80(140)	OE-HR 80	3
<p>● HR 86</p> 	G 312 /0	HR 86(160)	OE-HR 86	3










CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL

SPRING-EYED NEEDLES – NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
1/2-Kreis Rundkörper, fein = E(320) <i>1/2 circle, round bodied, thin = E(320)</i>				
● HR 18 F 	E 320/6	HR 18(60)	OE-HR 18F	6
● HR 20 F 	E 320/5	HR 20(60)	OE-HR 20F	6
● HR 23 F 	E 320/4	HR 23(70)	OE-HR 23F	6
● HR 54 F 	E 320/0	HR 54(100)	OE-HR 54F	3
● HR 63 F 	E 320/00	HR 63(110)	OE-HR 63F	3
● HR 72 F 	E 320/000	HR 72(110)	OE-HR 72F	3
1/2-Kreis Rundkörper, extra fein = Pb(332) <i>1/2 circle, round bodied, extra thin = Pb(332)</i>				
● HR 24 FF 	Pb 332/5	HR 24(60)	OE-HR 24FF	6
● HR 30 FF 	Pb 332/4	HR 30(60)	OE-HR 30FF	6
● HR 36 FF 	Pb 332/3	HR 36(60)	OE-HR 36FF	6
● HR 42 FF 	Pb 332/2	HR 42(70)	OE-HR 42FF	6












CHIRURGISCHE FEDERÖHRNADELN - UNSTERIL

SPRING-EYED NEEDLES - NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
1/2-Kreis Rundkörper, stark = GaR(314) <i>1/2 circle, round bodied, strong diameter = GaR(314)</i>				
● HR 23 S 	GaR 314 /8	HR 23(100)	OE-HR 23S	6
● HR 27 S 	GaR 314 /7	HR 27(100)	OE-HR 27S	6
● HR 32 S 	GaR 314 /6	HR 32(110)	OE-HR 32S	6
● HR 38 S 	GaR 314 /5	HR 38(130)	OE-HR 38S	6
● HR 46 S 	GaR 314 /4	HR 46(140)	OE-HR 46S	6
● HR 54 S 	GaR 314 /3	HR 54(140)	OE-HR 54S	3
● HR 60 S 	GaR 314 /2	HR 60(160)	OE-HR 60S	3
1/2-Kreis schneidende Rundkörper (Diamond Spitze) = G(412) <i>1/2 circle, round bodied, trocar point = G(412)</i>				
▼ HRT 20 	G 412 /15	HRT 20(70)	OE-HRT 20	6
▼ HRT 23 	G 412 /14	HRT 23(80)	OE-HRT 23	6











CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL

SPRING-EYED NEEDLES – NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
 ▼ HRT 26	G 412 /13	HRT 26(80)	OE-HRT 26	6
 ▼ HRT 30	G 412 /12	HRT 30(80)	OE-HRT 30	6
 ▼ HRT 34	G 412 /11	HRT 34(90)	OE-HRT 34	6
 ▼ HRT 38	G 412 /10	HRT 38(90)	OE-HRT 38	6
 ▼ HRT 42	G 412 /9	HRT 42(100)	OE-HRT 42	6
 ▼ HRT 46	G 412 /8	HRT 46(100)	OE-HRT 46	6
 ▼ HRT 50	G 412 /7	HRT 50(110)	OE-HRT 50	6
 ▼ HRT 54	G 412 /6	HRT 54(110)	OE-HRT 54	3
 ▼ HRT 58	G 412 /5	HRT 58(120)	OE-HRT 58	3
 ▼ HRT 63	G 412 /4	HRT 63(130)	OE-HRT 63	3
 ▼ HRT 68	G 412 /3	HRT 68(130)	OE-HRT 68	3










CHIRURGISCHE FEDERÖHRNADELN - UNSTERIL

SPRING-EYED NEEDLES - NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
Gerade schneidende Rundkörper (Diamond Spitze) <i>straight, round bodied, trocar point</i>				
▼ GRT 15 		GRT 15(80)	OE-GRT 15	20
1/2-Kreis außen schneidende, standard = GT (212) <i>1/2 circle, reverse cutting, standard = GT (212)</i>				
▼ HS 23 	GT 212 /14	HS 23(80)	OE-HS 23	6
▼ HS 26 	GT 212 /13	HS 26(80)	OE-HS 26	6
▼ HS 30 	GT 212 /12	HS 30(80)	OE-HS 30	6
▼ HS 34 	GT 212 /11	HS 34(90)	OE-HS 34	6
▼ HS 38 	GT 212 /10	HS 38(90)	OE-HS 38	6
▼ HS 42 	GT 212 /9	HS 42(100)	OE-HS 42	6
▼ HS 46 	GT 212 /8	HS 46(100)	OE-HS 46	6
▼ HS 50 	GT 212 /7	HS 50(110)	OE-HS 50	6
▼ HS 54 	GT 212 /6	HS 54(110)	OE-HS 54	3












CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL

SPRING-EYED NEEDLES – NON STERILE





Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
▼ HS 58 	GT 212 /5	HS 58(120)	OE-HS 58	3
▼ HS 63 	GT 212 /4	HS 63(130)	OE-HS 63	3
▼ HS 68 	GT 212 /3	HS 68(130)	OE-HS 68	3
▼ HS 73 	GT 212 /2	HS 73(140)	OE-HS 73	3
<p>1/2-Kreis außen schneidende, stark = GaT (214) <i>1/2 circle, reverse cutting, strong diameter = GaT(214)</i></p>				
▼ HS 18 S 	GaT 214 /9	HS 18(90)	OE-HS 18S	6
▼ HS 23 S 	GaT 214 /8	HS 23(100)	OE-HS 23S	6
▼ HS 27 S 	GaT 214 /7	HS 27(100)	OE-HS 27S	6
▼ HS 32 S 	GaT 214 /6	HS 32(110)	OE-HS 32S	6
▼ HS 38 S 	GaT 214 /5	HS 38(130)	OE-HS 38S	6

CHIRURGISCHE FEDERÖHRNADELN - UNSTERIL

SPRING-EYED NEEDLES - NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
▼ HS 46 S 	GaT 214 /4	HS 46(140)	OE-HS 46S	6
▼ HS 54 S 	GaT 214 /3	HS 54(140)	OE-HS 54S	3
3/8-Kreis außen schneidende, standard = BT (204) <i>3/8 circle, reverse cutting, standard = BT (204)</i>				
▼ DS 18 	BT 204 /16	DS 18(70)	OE-DS 18	6
▼ DS 20 	BT 204 /15	DS 20(70)	OE-DS 20	6
▼ DS 23 	BT 204 /14	DS 23(80)	OE-DS 23	6
▼ DS 26 	BT 204 /13	DS 26(80)	OE-DS 26	6
▼ DS 30 	BT 204 /12	DS 30(80)	OE-DS 30	6
▼ DS 34 	BT 204 /11	DS 34(90)	OE-DS 34	6
▼ DS 38 	BT 204 /10	DS 38(90)	OE-DS 38	6
▼ DS 42 	BT 204 /9	DS 42(100)	OE-DS 42	6
▼ DS 46 	BT 204 /8	DS 46(100)	OE-DS 46	3

CHIRURGISCHE FEDERÖHRNADELN – UNSTERIL SPRING-EYED NEEDLES – NON STERILE

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Inhalt/Blister blister content
▼ DS 50 	BT 204 /7	DS 50(110)	OE-DS 50	3
▼ DS 54 	BT 204 /6	DS 54(110)	OE-DS 54	3
▼ DS 70 	BT 204 /3	DS 70(130)	OE-DS 70	3
Gerade schneidende, standard (Keith) <i>straight cutting, standard (Keith)</i>				
▼ GS 20 		GS 20(80)	OE-GS 20	20